

### **2,3-Dichlorobiphenyl (PCB 5) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007066  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter                | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|----------------------------------|--------------|---------------------|-------|--------|--------|
|                  | 2,3-Dichlorobiphenyl (P<br>CB 5) | [16605-91-7] | 100                 | µg/ml |        |        |