

2,4-Dichlorobiphenyl (PCB 7) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007073
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,4-Dichlorobiphenyl (P CB 7) | [33284-50-3] | 500 | µg/ml | | |