

### **3,4-Dichlorobiphenyl (PCB 12) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007087  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4-Dichlorobiphenyl (P CB 12)	[2974-92-7]	100	µg/ml		