

3,4-Dichlorobiphenyl (PCB 12) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007088
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4-Dichlorobiphenyl (P CB 12)	[2974-92-7]	500	µg/ml		