

### **3,4'-Dichlorobiphenyl (PCB 13) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007090  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4'-Dichlorobiphenyl ( PCB 13)	[2974-90-5]	100	µg/ml		