

### **3,4'-Dichlorobiphenyl (PCB 13) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007091  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter                  | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|------------------------------------|-------------|---------------------|-------|--------|--------|
|                  | 3,4'-Dichlorobiphenyl (<br>PCB 13) | [2974-90-5] | 500                 | µg/ml |        |        |