

### **4,4'-Dichlorobiphenyl (PCB 15) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007097  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4,4'-Dichlorobiphenyl ( PCB 15)	[2050-68-2]	500	µg/ml		