

2,3,4'-Trichlorobiphenyl (PCB 22) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007105
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,3,4'-Trichlorobipheny I (PCB 22) | [38444-85-8] | 100 | µg/ml | | |