

2,3',4-Trichlorobiphenyl (PCB 25) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007114
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',4-Trichlorobiphenyl I (PCB 25)	[55712-37-3]	100	µg/ml		