

## **2,3',6-Trichlorobiphenyl (PCB 27) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007120  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter                      | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
|                  | 2,3',6-Trichlorobiphenyl<br>I (PCB 27) | [38444-76-7] | 100                 | µg/ml |        |        |