

2',3,4-Trichlorobiphenyl (PCB 33) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007138
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,4-Trichlorobiphenyl I (PCB 33)	[38444-86-9]	100	µg/ml		