

3,4,4'-Trichlorobiphenyl (PCB 37) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007150
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4,4'-Trichlorobiphenyl I (PCB 37)	[38444-90-5]	100	µg/ml		