

## **2,2',5,5'-Tetrachlorobiphenyl (PCB 52) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007159  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',5,5'-Tetrachlorobi phenyl (PCB 52)	[35693-99-3]	100	µg/ml		