

2,3,3',4'-Tetrachlorobiphenyl (PCB 56) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007165
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
| | 2,3,3',4'-Tetrachlorobi phenyl (PCB 56) | [41464-43-1] | 100 | µg/ml | | |