

## **2,3,3',5'-Tetrachlorobiphenyl (PCB 58) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007172  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',5'-Tetrachlorobi phenyl (PCB 58)	[41464-49-7]	500	µg/ml		