

## **2,3,4',5-Tetrachlorobiphenyl (PCB 63)**

Art. ID                      LM24-OP-4007179  
Unit                         10 mg  
Deliverydetails            De Minimis UN3432 II E2

### Description

Shelf-Life: 36 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,4',5-Tetrachlorobiphenyl (PCB 63)	[74472-34-7]				