

## **Low alloyed white cast iron - chips**

Art. ID                      ECRM-F 488-2  
Unit                         100 g  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,956	%		
	Silicon (Si)	[7440-21-3]	0,374	%		
	Manganese (Mn)	[7439-96-5]	0,201	%		
	Phosphorus (P)	[7723-14-0]	0,0111	%		
	Sulfur (S)	[7704-34-9]	0,1173	%		
	Chromium (Cr)	[7440-47-3]	0,303	%		
	Nickel (Ni)	[7440-02-0]	0,1247	%		
	Copper (Cu)	[7440-50-8]	0,0256	%		
	Nitrogen (N)	[7727-37-9]	0,0052	%		
	Tin (Sn)	[7440-31-5]	0,0013	%		
	Titanium (Ti)	[7440-32-6]	0,0636	%		
	Vanadium (V)	[7440-62-2]	0,0545	%		