

## **2,3',4',6-Tetrachlorobiphenyl (PCB 71) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007199  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',4',6-Tetrachlorobi phenyl (PCB 71)	[41464-46-4]	100	µg/ml		