

## **2,3',5',6-Tetrachlorobiphenyl (PCB 73) 500 µg/ml in Iso-octane**

Art. ID                      LM24-OP-4007206  
Unit                         1 mL  
Deliverydetails            De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',5',6-Tetrachlorobi phenyl (PCB 73)	[74338-23-1]	500	µg/ml		