

2,4,4',5-Tetrachlorobiphenyl (PCB 74) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007209
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|-------|--------|--------|
| | 2,4,4',5-Tetrachlorobiphenyl (PCB 74) | [32690-93-0] | 500 | µg/ml | | |