

**2,4,4',6-Tetrachlorobiphenyl (PCB 75) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007211  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,4',6-Tetrachlorobiphenyl (PCB 75)	[32598-12-2]	100	µg/ml		