

2',3,4,5-Tetrachlorobiphenyl (PCB 76) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007215
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,4,5-Tetrachlorobiphenyl (PCB 76)	[70362-48-0]	500	µg/ml		