

3,3',4,4'-Tetrachlorobiphenyl (PCB 77)

Art. ID LM24-OP-4007216
Unit 20 mg
Deliverydetails De Minimis UN3432 II E2

Description

Shelf-Life: 36 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|------|--------|--------|
| | 3,3',4,4'-Tetrachlorobi phenyl (PCB 77) | [32598-13-3] | | | | |