

3,3',4,4'-Tetrachlorobiphenyl (PCB 77) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007218
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4'-Tetrachlorobi phenyl (PCB 77)	[32598-13-3]	500	µg/ml		