

### **3,3',4,5-Tetrachlorobiphenyl (PCB 78) 500 µg/ml in Iso-octane**

Art. ID                      LM24-OP-4007221  
Unit                         1 mL  
Deliverydetails            De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,5-Tetrachlorobiphenyl (PCB 78)	[70362-49-1]	500	µg/ml		