

3,3',4,5'-Tetrachlorobiphenyl (PCB 79) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007223
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,5'-Tetrachlorobi phenyl (PCB 79)	[41464-48-6]	100	µg/ml		