

3,4,4',5-Tetrachlorobiphenyl (PCB 81) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007230
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4,4',5-Tetrachlorobiphenyl (PCB 81)	[70362-50-4]	100	µg/ml		