

### **3,4,4',5-Tetrachlorobiphenyl (PCB 81) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007231  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,4,4',5-Tetrachlorobiphenyl (PCB 81)	[70362-50-4]	500	µg/ml		