

## **2,2',4,5,5'-Pentachlorobiphenyl (PCB 101) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007234  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',4,5,5'-Pentachloro biphenyl (PCB 101)	[37680-73-2]	500	µg/ml		