

## **2,3,3',4,5'-Pentachlorobiphenyl (PCB 108) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007242  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',4,5'-Pentachloro biphenyl (PCB 108)	[70362-41-3]	100	µg/ml		