

2,3,3',5,5'-Pentachlorobiphenyl (PCB 111) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007249
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',5,5'-Pentachloro biphenyl (PCB 111)	[39635-32-0]	500	µg/ml		