

## **2,3',4,4',5-Pentachlorobiphenyl (PCB 118) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007252  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter                             | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|---|--------------|---------------------|-------|--------|--------|
|                  | 2,3',4,4',5-Pentachloro<br>biphenyl (PCB 118) | [31508-00-6] | 100                 | µg/ml |        |        |