

## **2,3',4,4',6-Pentachlorobiphenyl (PCB 119) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007256  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',4,4',6-Pentachloro biphenyl (PCB 119)	[56558-17-9]	500	µg/ml		