

**2',3,4,4',5-Pentachlorobiphenyl (PCB 123) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007268  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2',3,4,4',5-Pentachloro biphenyl (PCB 123)	[65510-44-3]	500	µg/ml		