

3,3',4,5,5'-Pentachlorobiphenyl (PCB 127) 100 µg/ml in Iso-octane

Art. ID LM24-OP-4007280
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,5,5'-Pentachloro biphenyl (PCB 127)	[39635-33-1]	100	µg/ml		