

2,3,3',4,4',5'-Hexachlorobiphenyl (PCB 157) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007287
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3,3',4,4',5'-Hexachlorobiphenyl (PCB 157)	[69782-90-7]	500	µg/ml		