

## **2,3',4,4',5',6-Hexachlorobiphenyl (PCB 168) 500 µg/ml in Iso-octane**

Art. ID LM24-OP-4007297  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,3',4,4',5',6-Hexachlorobiphenyl (PCB 168)	[59291-65-5]	500	µg/ml		