

### **3,3',4,4',5,5'-Hexachlorobiphenyl (PCB 169) 100 µg/ml in Iso-octane**

Art. ID LM24-OP-4007300  
Unit 1 mL  
Deliverydetails De Minimis UN1262 II E2

#### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4',5,5'-Hexachlorobiphenyl (PCB 169)	[32774-16-6]	100	µg/ml		