

3,3',4,4',5,5'-Hexachlorobiphenyl (PCB 169) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007301
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',4,4',5,5'-Hexachlorobiphenyl (PCB 169)	[32774-16-6]	500	µg/ml		