

2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180) 500 µg/ml in Iso-octane

Art. ID LM24-OP-4007304
Unit 1 mL
Deliverydetails De Minimis UN1262 II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,2',3,4,4',5,5'-Heptachlorobiphenyl (PCB 180)	[35065-29-3]	500	µg/ml		