

HighPurS-CERT - Potassium Nitrate (KNO₃) - min. 99.999%, ISO 17034 and ISO/IEC 17025 certified

Art. ID LM24-OP-4007350
Unit 50 g
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1486 III E1

Description

Shelf-life: 36 months

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium nitrate	[7757-79-1]				