

HighPurS-CERT - Manganese(II) nitrate (Mn(NO₃)₂*H₂O) - min. 99.99+%, ISO 17034 and ISO/IEC 17025 certified

Art. ID LM24-OP-4007352
Unit 100 g
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) to be checked

Description

Shelf-Life: 36 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Manganese(II)nitrate					