

## **DNPH-Mix - 15 components, 100 µg/ml each of in Acetonitrile**

Art. ID LM24-OP-4007378  
Unit 1 mL  
Deliverydetails De Minimis UN1648 II E2

### Description

Shelf-Life: 6 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acetaldehyde-2,4-dinitrophenylhydrazone	[1019-57-4]	100	µg/mL		
	Acetone 2,4-Dinitrophenylhydrazone	[1567-89-1]	100	µg/mL		
	Acrolein 2,4-dinitrophenylhydrazone	[888-54-0]	100	µg/mL		
	Benzaldehyde-DNPH	[1157-84-2]	100	µg/mL		
	Butyraldehyde 2,4-Dinitrophenylhydrazone	[1527-98-6]	100	µg/mL		
	Crotonaldehyde-2,4-dinitrophenylhydrazone	[1527-96-4]	100	µg/mL		
	Cyclohexanon-DNPH	[1589-62-4]	100	µg/mL		
	Formaldehyde-2,4-dinitrophenylhydrazone (DNPH)	[1081-15-8]	100	µg/mL		
	Isovaleraldehyde-DNPH	[2256-01-1]	100	µg/mL		
	m-Tolualdehyde-DNPH	[2880-05-9]	100	µg/mL		
	o-Tolualdehyde-DNPH	[1773-44-0]	100	µg/mL		
	p-Tolualdehyde-DNPH	[2571-00-8]	100	µg/mL		
	Valeraldehyde-DNPH	[2057-84-3]	100	µg/mL		
	Propionaldehyde 2,4-Dinitrophenylhydrazone	[725-00-8]	100	µg/mL		