

### **ICH Q3D parenteral, Standard 1 - 10 components in 12% HNO3**

Art. ID LM24-OP-4007419  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) to be checked

#### **Description**

Shelf-Life: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	10	mg/l		
	Arsenic (As)	[7440-38-2]	15	mg/l		
	Cadmium (Cd)	[7440-43-9]	2	mg/l		
	Cobalt (Co)	[7440-48-4]	5	mg/l		
	Mercury (Hg)	[7439-97-6]	3	mg/l		
	Nickel (Ni)	[7440-02-0]	20	mg/l		
	Lead (Pb)	[7439-92-1]	5	mg/l		
	Selenium (Se)	[7782-49-2]	80	mg/l		
	Thallium (Tl)	[7440-28-0]	8	mg/l		
	Vanadium (V)	[7440-62-2]	10	mg/l		