

### **ICH Q3D parenteral, Standard 3 - 7 components in 5% HNO<sub>3</sub>/0,2% HF**

Art. ID LM24-OP-4007420  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

#### Description

Shelf-Life: 12 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	70	mg/l		
	Chromium (Cr)	[7440-47-3]	110	mg/l		
	Copper (Cu)	[7440-50-8]	30	mg/l		
	Lithium (Li)	[7439-93-2]	25	mg/l		
	Molybdenum (Mo)	[7439-98-7]	150	mg/l		
	Antimony (Sb)	[7440-36-0]	9	mg/l		
	Tin (Sn)	[7440-31-5]	60	mg/l		