

USP 232/ ICH Q3D Elemental Impurities Standard 2-1: 6 components in 2% HNO₃

Art. ID LM24-OP-4007421
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cobalt (Co)	[7440-48-4]	50	µg/ml		
	Vanadium (V)	[7440-62-2]	100	µg/ml		
	Nickel (Ni)	[7440-02-0]	200	µg/ml		
	Thallium (Tl)	[7440-28-0]	8	µg/ml		
	Selenium (Se)	[7782-49-2]	150	µg/ml		
	Silver (Ag)	[7440-22-4]	150	µg/ml		