

USP 232/ ICH Q3D Elemental Impurities Standard 3 - 7 components in 5% HNO3/0,1% HF

Art. ID LM24-OP-4007422
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Lithium (Li)	[7439-93-2]	550	µg/ml		
	Antimony (Sb)	[7440-36-0]	1200	µg/ml		
	Barium (Ba)	[7440-39-3]	1400	µg/ml		
	Molybdenum (Mo)	[7439-98-7]	3000	µg/ml		
	Copper (Cu)	[7440-50-8]	3000	µg/ml		
	Tin (Sn)	[7440-31-5]	6000	µg/ml		
	Chromium (Cr)	[7440-47-3]	11000	µg/ml		