

USP 232/ ICH Q3D Elemental Impurities Standard 3 - 7 components in 5% HNO3/0,1% HF

Art. ID LM24-OP-4007422
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2922 (6.1) II E2

Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Lithium (Li) | [7439-93-2] | 550 | µg/ml | | |
| | Antimony (Sb) | [7440-36-0] | 1200 | µg/ml | | |
| | Barium (Ba) | [7440-39-3] | 1400 | µg/ml | | |
| | Molybdenum (Mo) | [7439-98-7] | 3000 | µg/ml | | |
| | Copper (Cu) | [7440-50-8] | 3000 | µg/ml | | |
| | Tin (Sn) | [7440-31-5] | 6000 | µg/ml | | |
| | Chromium (Cr) | [7440-47-3] | 11000 | µg/ml | | |