

Zinc (Zn) 10000 mg/L in 2-5% HNO₃, Single-Element Standard for AAS (Flame)

Art. ID LM24-OP-4007460
Unit 20 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Zinc (Zn) | [7440-66-6] | 10000 | mg/L | | |