

## **2-Heptanone**

Art. ID                    LM24-OP-4007507  
Unit                      1 mL  
Deliverydetails        De Minimis UN1110 III E1

### Description

Shelf-Life: 24 months /// Produced under ISO 17034 and ISO 17025

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | 2-Heptanone       | [110-43-0] |                     |      |        |        |